



- **Portable thromboelastograph** with **immediate** response
- The only instrument for thromboelastography that can be used directly **near the animal**
- The thromboelastography test cannot be outsourced and therefore this is the best solution for performing the test in the veterinary clinic
- **Minimum sample volume: 230 µl** using the capillary integrated in the cartridge
- **The only truly system** based on **whole blood samples** with **no chemical activators**
- **Instant first response** in *hypercoagulant* patients
- **Fast first response times: 10 min** for initial assessment of the coagulation process (CT, CFT phases, alpha angle) for *healty/ hypogoagulant* patients
- Power supply via **fixed electric network** or via **integrated battery**
- Electronic and external controls available
- **No maintenance** required



The VCM system provides the following parameters:

CT: clot time (sec): time from the start of the test until a height of 1% above the baseline is reached

CFT: clot formation time (sec): time from 1% to 10% height on baseline

Alpha (°): is the angle between the horizontal axis (times) and the tangent to the curve at 1% amplitude. Describes the kinetics of coagulation

A10, A20 (VCM units): amplitude at 10 and 20 min. They represent the robustness of the clot 10 and 20 minutes after its formation

MCF (VCM unit): maximum clot formation: it is the measure of clot strength and therefore of its quality. It corresponds to the maximum amplitude of the clot before the start of the fibrinolytic phase

LI30, LI45 (%): is the clot size at 30 and 45 min as a percentage of the MCF

COMPLETE REPORT

OF ALL QUANTITATIVE PARAMETERS (CT, CFT, ALPHA ANGLE, MCF, LI30, LI45)

PORTABLE SYSTEM

CAN BE USED DIRECTLY ON THE ANIMAL

BATTERY ON-BOARD

OR MAINS POWER SUPPLY

FRESH WHOLE BLOOD

SAMPLES WITH NO CHEMICAL ACTIVATORS

Entegrión VCM VET


PORTABLE SYSTEM FOR THE FAST AND COMPLETE
ASSESSMENT OF THE GLOBAL HEMOSTATIC CONDITION



Some examples of layouts obtained with the **VCM Vet** system:

SISTEMA > RESULT


CT	7.9m	[5.2-8.7]
CFT	2.8m	[1.7-3.9]
Alpha	53°	[46-65]
A10	26	[19-33]
A20	34	[25-41]
MCF	37	[27-44]
LI30	100%	[97-100]
LI45	99%	[89-100]



Healthy profile

SISTEMA > RESULT


CT	6.2m	[5.2-8.7]
CFT	1.6m	[1.7-3.9]
Alpha	68°	[46-65]
A10	35	[19-33]
A20	42	[25-41]
MCF	45	[27-44]
LI30	100%	[97-100]
LI45	100%	[89-100]



From a hypercoagulant animal

SISTEMA > RESULT


CT	16.0m	[5.2-8.7]
CFT	9.4m	[1.7-3.9]
Alpha	25°	[46-65]
A10	11	[19-33]
A20	16	[25-41]
MCF	19	[27-44]
LI30	100%	[97-100]
LI45		



From a hypocoagulant animal

SYSTEM > RESULT

CT	28.5m	[5.2-8.7]
CFT	22.6m	[1.7-3.9]
Alpha	9°	[46-65]
A10	4	[19-33]
A20	9	[25-41]
MCF	15	[27-44]
LI30	100%	[97-100]
LI45		



Effect of UFH administration
(unfractionated heparin)

Code Description

Analyzer:

9001PCM VCM Vet System

Test:

3105PCM VCM Vet Test Cartridge Pack

Accessories:

3103PCM VCM Vet Heater plate Pack Insert ML

9011PCM VCM Vet Level 1 Quality Control kit

9012PCM VCM Vet Level 2 Quality Control kit

Eos srl

Via Fossona 7/A, 35030 Cervarese S. Croce (PD), Italia
tel. +39 049 991 5554; fax +39 049 991 9469
info@eosbio.com; www.eosbio.com